


Product Name: Recombinant Human GLI1 Protein
Catalog #: PEH4775



Summary

Name	GLI1_HUMAN
Purity	Greater than 85% as determined by reducing SDS-PAGE
Endotoxin level	Please contact with the lab for this information
Tag	C-His
Accession #	P08151
Host	E.coli
Species	Human
Molecular Mass	27 kDa
Buffer	Phosphate buffered saline
Shipping	The product is shipped on dry ice/polar packs. Upon receipt, store it immediately at the temperature listed below.
Stability&Storage	Store at $\leq -70^{\circ}\text{C}$, stable for 6 months after receipt. Store at $\leq -70^{\circ}\text{C}$, stable for 3 months under sterile conditions after opening. Please minimize freeze-thaw cycles.
Sequence	AFSNASDRAK HQRNTHSNEK PYVCKLPGCT KRYTDPSSLR KHVKTVHGPD AHVTKRHRGD GPLPRAPRIS TVEPKREREG GPIREESRLT VPEGAMKPQP SPGAQSSCSS DHSPAGSAAN TDSGVEMTGN AGGSTEDLSS LDEGPCIAGT GLSTLRRLEN LRLDQLHQLR PIGTRGLKLP SLSHTGTTVS RRVGPPVSLE RRSSSSSSIS SAYTVSRRSS LASPFPPGSP

image



Background

Alternative Names Zinc finger protein GLI1 (Glioma-associated oncogene) (Oncogene GLI)

Note

Product Name: Recombinant Human GLI1 Protein
Catalog #: PEH4775



For Research Use Only , Not for Diagnostic Use