


Product Name: Recombinant Human USP17L7 Protein
Catalog #: PEH4777



Summary

Name	U17L7_HUMAN
Purity	Greater than 85% as determined by reducing SDS-PAGE
Endotoxin level	Please contact with the lab for this information
Tag	C-His
Accession #	P0C7H9
Host	E.coli
Species	Human
Molecular Mass	15 kDa
Buffer	Lyophilized from Phosphate buffered saline
Shipping	The product is shipped on dry ice/polar packs. Upon receipt, store it immediately at the temperature listed below.
Stability&Storage	Store at $\leq -70^{\circ}\text{C}$, stable for 12 months after receipt. Store at $\leq -70^{\circ}\text{C}$, stable for 3 months under sterile conditions after opening. Please minimize freeze-thaw cycles.
Sequence	RPATQGELKR DHPCLQVPEL DEHLVERATQ ESTLDHWKFP QEQNKTKPEF NVRKVEGTLP PNVLVIHQSK YKCGMKNHHP EQQSSLLNLS STKPTDQESM NTGTLASLQG STRRSKGNNK

image



Background

Alternative Names Inactive ubiquitin carboxyl-terminal hydrolase 17-like protein 7

Note

For Research Use Only , Not for Diagnostic Use