


Product Name: Recombinant Human ZNF286B Protein
Catalog #: PEH4782



Summary

Name	Z286B_HUMAN
Purity	Greater than 85% as determined by reducing SDS-PAGE
Endotoxin level	Please contact with the lab for this information
Tag	C-His
Accession #	P0CG31
Host	E.coli
Species	Human
Molecular Mass	17 kDa
Buffer	Lyophilized from Phosphate buffered saline
Shipping	The product is shipped on dry ice/polar packs. Upon receipt, store it immediately at the temperature listed below.
Stability&Storage	Store at $\leq -70^{\circ}\text{C}$, stable for 12 months after receipt. Store at $\leq -70^{\circ}\text{C}$, stable for 3 months under sterile conditions after opening. Please minimize freeze-thaw cycles.
Sequence	PSDVETRPQS KDSTSVQDFS KAESCKVAII DRLTRNSVYD SNLEAALECE NWLEKQQGNQ ERHLREMFTH MNSLSEETHD EHDVYWKSFN QKSVLITEDR VPKGSYAFHT LEKSLKQKSN

image



Background

Alternative Names Putative zinc finger protein 286B

Note

For Research Use Only , Not for Diagnostic Use