

Product Name: Recombinant Human ERVK-10 Protein
Catalog #: PEH4786



Summary

Name	VPK10_HUMAN
Purity	Greater than 85% as determined by reducing SDS-PAGE
Endotoxin level	Please contact with the lab for this information
Tag	C-His
Accession #	P10265
Host	E.coli
Species	Human
Molecular Mass	18 kDa
Buffer	Phosphate buffered saline
Shipping	The product is shipped on dry ice/polar packs. Upon receipt, store it immediately at the temperature listed below.
Stability&Storage	Store at $\leq -70^{\circ}\text{C}$, stable for 6 months after receipt. Store at $\leq -70^{\circ}\text{C}$, stable for 3 months under sterile conditions after opening. Please minimize freeze-thaw cycles.
Sequence	VCKAIQ GKQ FEGLVDTGAD VSIILNQWP KNWPKQKAVT GLVGIGTASE VYQSMEILHC LGPDNQESTV QPMITSIPLN LWGRDLLQQW GAEITMPAPL YSPTSQKIMT KMGYIPGKGL GKNEDGIKVP VEA KINQERE

image



Background

Alternative Names	Endogenous retrovirus group K member 10 Pro protein (HERV-K10 Pro protein) (HERV-K107 Pro protein) (HERV-K_5q33.3 provirus ancestral Pro protein) (EC 3.4.23.50) (Protease) (Proteinase) (PR)
--------------------------	---

Note

For Research Use Only , Not for Diagnostic Use

Product Name: Recombinant Human ERVK-10 Protein
Catalog #: PEH4786

