


Product Name: Recombinant Human MYCLP1 Protein
Catalog #: PEH4789



Summary

Name	MYCP1_HUMAN
Purity	Greater than 85% as determined by reducing SDS-PAGE
Endotoxin level	Please contact with the lab for this information
Tag	C-His
Accession #	P12525
Host	E.coli
Species	Human
Molecular Mass	16 kDa
Buffer	Phosphate buffered saline
Shipping	The product is shipped on dry ice/polar packs. Upon receipt, store it immediately at the temperature listed below.
Stability&Storage	Store at $\leq -70^{\circ}\text{C}$, stable for 6 months after receipt. Store at $\leq -70^{\circ}\text{C}$, stable for 3 months under sterile conditions after opening. Please minimize freeze-thaw cycles.
Sequence	FATPDYTPEL EAGNLAPIFP CLLGEPKIQA CSRSESPSDS EGEEIDVTVK KRQSLSTRKP VIIAVRADLL DPRMNLFHIS IHQQQHNYAA PFPPESCFQE GAPKRMPPKE ALEREAPGGK DDKEDDEEIVS

image



Background

Alternative Names	Putative myc-like protein MYCLP1 (Protein L-Myc-2) (V-myc myelocytomatosis viral oncogene homolog pseudogene 1)
--------------------------	---

Note

For Research Use Only , Not for Diagnostic Use