



Product Name: Recombinant Human ZNF35 Protein
Catalog #: PEH4790



Summary

Name	ZNF35_HUMAN
Purity	Greater than 85% as determined by reducing SDS-PAGE
Endotoxin level	Please contact with the lab for this information
Tag	C-His
Accession #	P13682
Host	E.coli
Species	Human
Molecular Mass	23 kDa
Buffer	Phosphate buffered saline
Shipping	The product is shipped on dry ice/polar packs. Upon receipt, store it immediately at the temperature listed below.
Stability&Storage	Store at $\leq -70^{\circ}\text{C}$, stable for 6 months after receipt. Store at $\leq -70^{\circ}\text{C}$, stable for 3 months under sterile conditions after opening. Please minimize freeze-thaw cycles.
Sequence	KEEEEEENFP GQASSQQVHS ENIKVWAPVQ GLQTGLDGSE EEEKGQNISW DMAVVLKATQ EAPAASTLGS YSLPGTLAKS EILETHGTMN FLGAETKNLQ LLVPKTEICE EAEKPLIISE RIQKADPQGP ELGEACEKGN MLKRQRIKRE KKDFRQVIVN DCHLPESFKE EENQKCKKSG GKYSLSNGAV

image



Background

Alternative Names Zinc finger protein 35 (Zinc finger protein HF.10)

Note

For Research Use Only , Not for Diagnostic Use

Product Name: Recombinant Human ZNF35 Protein
Catalog #: PEH4790

