


Product Name: Recombinant Human ZNF44 Protein
Catalog #: PEH4793



Summary

Name	ZNF44_HUMAN
Purity	Greater than 85% as determined by reducing SDS-PAGE
Endotoxin level	Please contact with the lab for this information
Tag	C-His
Accession #	P15621
Host	E.coli
Species	Human
Molecular Mass	15 kDa
Buffer	Phosphate buffered saline
Shipping	The product is shipped on dry ice/polar packs. Upon receipt, store it immediately at the temperature listed below.
Stability&Storage	Store at $\leq -70^{\circ}\text{C}$, stable for 6 months after receipt. Store at $\leq -70^{\circ}\text{C}$, stable for 3 months under sterile conditions after opening. Please minimize freeze-thaw cycles.
Sequence	GDGPHKCKIC GKGDFPGSA RIHEGHTLE KPYECKQCGK LLSHRSSFRR HMMAHTGDGP HKCTVCGKAF DSPSVFQRHE RTHTGEKPYE CKQCGKAFRT SSSLRKHETT HTGEQPYKCK

image



Background

Alternative Names Zinc finger protein 44 (Gonadotropin-inducible ovary transcription repressor 2) (GIOT-2) (Zinc finger protein 55) (Zinc finger protein 58) (Zinc finger protein KOX7)

Note

For Research Use Only , Not for Diagnostic Use

Product Name: Recombinant Human ZNF44 Protein
Catalog #: PEH4793

