


Product Name: Recombinant Human ZNF24 Protein
Catalog #: PEH4796



Summary

Name	ZNF24_HUMAN
Purity	Greater than 85% as determined by reducing SDS-PAGE
Endotoxin level	Please contact with the lab for this information
Tag	C-His
Accession #	P17028
Host	E.coli
Species	Human
Molecular Mass	27 kDa
Buffer	Phosphate buffered saline
Shipping	The product is shipped on dry ice/polar packs. Upon receipt, store it immediately at the temperature listed below.
Stability&Storage	Store at $\leq -70^{\circ}\text{C}$, stable for 6 months after receipt. Store at $\leq -70^{\circ}\text{C}$, stable for 3 months under sterile conditions after opening. Please minimize freeze-thaw cycles.
Sequence	ILIIPTDEE EKILRVKLEE DPDGEEGSSI PWNHLPDPEI FRQRFQFGY QDSPGPRAEV SQLRELCRLW LRPETHTEQ ILELVVLEQF VAILPKELQT WVRDHPENG EEAFTVLEDL ESELDGPGQP VSLRRRKREV LVEDMVSQEE AQGLPSSELD AVENQLKQAS WELHSLRHCD DDGRTENGAL APKQELPSAL ESHEVPGTLN MGVPQIFKYG

image



Background

Alternative Names	Zinc finger protein 24 (Retinoic acid suppression protein A) (RSG-A) (Zinc finger and SCAN domain-containing protein 3) (Zinc finger protein 191) (Zinc finger protein KOX17)
--------------------------	---

Product Name: Recombinant Human ZNF24 Protein
Catalog #: PEH4796



Note

For Research Use Only , Not for Diagnostic Use