


Product Name: Recombinant Human ZKSCAN1 Protein
Catalog #: PEH4797



Summary

Name	ZKSC1_HUMAN
Purity	Greater than 85% as determined by reducing SDS-PAGE
Endotoxin level	Please contact with the lab for this information
Tag	C-His
Accession #	P17029
Host	E.coli
Species	Human
Molecular Mass	21 kDa
Buffer	Phosphate buffered saline
Shipping	The product is shipped on dry ice/polar packs. Upon receipt, store it immediately at the temperature listed below.
Stability&Storage	Store at $\leq -70^{\circ}\text{C}$, stable for 6 months after receipt. Store at $\leq -70^{\circ}\text{C}$, stable for 3 months under sterile conditions after opening. Please minimize freeze-thaw cycles.
Sequence	GLSPQAAQEK DGIVIVKVEE EDEEDHMWGQ DSTLQDTPPP DPEIFRQRFR RFCYQNTFGP REALSRLKEL CHQWLRPEIN TKEQILELLV LEQFLSILPK ELQVWLQEYR PDSGEEAVTL LEDLELDLSG QQVPGQVHGP EMLARGMVPL DPVQESSFD

image



Background

Alternative Names	Zinc finger protein with KRAB and SCAN domains 1 (Zinc finger protein 139) (Zinc finger protein 36) (Zinc finger protein KOX18)
--------------------------	---

Note

For Research Use Only , Not for Diagnostic Use

Product Name: Recombinant Human ZKSCAN1 Protein
Catalog #: PEH4797

