


Product Name: Recombinant Human ZNF30 Protein
Catalog #: PEH4799



Summary

Name	ZNF30_HUMAN
Purity	Greater than 85% as determined by reducing SDS-PAGE
Endotoxin level	Please contact with the lab for this information
Tag	C-His
Accession #	P17039
Host	E.coli
Species	Human
Molecular Mass	18 kDa
Buffer	Phosphate buffered saline
Shipping	The product is shipped on dry ice/polar packs. Upon receipt, store it immediately at the temperature listed below.
Stability&Storage	Store at $\leq -70^{\circ}\text{C}$, stable for 6 months after receipt. Store at $\leq -70^{\circ}\text{C}$, stable for 3 months under sterile conditions after opening. Please minimize freeze-thaw cycles.
Sequence	DDTIGCKEMP TSENCPSFAL HQKISRQKPR ECQEYGKTLC QDSKPVQHER IHSEKPNRC KECGKNFSNG HQLTIHQRLH VGEKPYKYEK CGKAFISGSA FVKHGRIHTG EKPLKCKQCG KTISGSYQLT VHKSHTGKK

image



Background

Alternative Names Zinc finger protein 30 (Zinc finger protein KOX28)

Note

For Research Use Only , Not for Diagnostic Use