


Product Name: Recombinant Human ZSCAN20 Protein
Catalog #: PEH4800



Summary

Name	ZSC20_HUMAN
Purity	Greater than 85% as determined by reducing SDS-PAGE
Endotoxin level	Please contact with the lab for this information
Tag	C-His
Accession #	P17040
Host	E.coli
Species	Human
Molecular Mass	23 kDa
Buffer	Phosphate buffered saline
Shipping	The product is shipped on dry ice/polar packs. Upon receipt, store it immediately at the temperature listed below.
Stability&Storage	Store at $\leq -70^{\circ}\text{C}$, stable for 6 months after receipt. Store at $\leq -70^{\circ}\text{C}$, stable for 3 months under sterile conditions after opening. Please minimize freeze-thaw cycles.
Sequence	EAWGEVANED AVKPSTLCPK APDMGFEMRH EDEDQISEQD IFEGLPGALS KCPTEAVCQP LDWGEDSENE NEDEGQWGNP SQEQWQESS EEDLEKLIDH QGLYLAEKPY KCDTCMKFS RSHFIAHQR IHTGEKPYKC LECGKNFSDR SNLNTHQRIH TGEKPYKCLE CGKSFS DHSN

image



Background

Alternative Names Zinc finger and SCAN domain-containing protein 20 (Zinc finger protein 31) (Zinc finger protein 360) (Zinc finger protein KOX29)

Note

For Research Use Only , Not for Diagnostic Use

Product Name: Recombinant Human ZSCAN20 Protein
Catalog #: PEH4800

