


Product Name: Recombinant Human EGR1 Protein
Catalog #: PEH4802



Summary

Name	EGR1_HUMAN
Purity	Greater than 85% as determined by reducing SDS-PAGE
Endotoxin level	Please contact with the lab for this information
Tag	C-His
Accession #	P18146
Host	E.coli
Species	Human
Molecular Mass	10 kDa
Buffer	Phosphate buffered saline
Shipping	The product is shipped on dry ice/polar packs. Upon receipt, store it immediately at the temperature listed below.
Stability&Storage	Store at $\leq -70^{\circ}\text{C}$, stable for 6 months after receipt. Store at $\leq -70^{\circ}\text{C}$, stable for 3 months under sterile conditions after opening. Please minimize freeze-thaw cycles.
Sequence	DHLTTHIRTH TGEKPFACDI CGRKFARSDE RKRHTKIHLR QKDKKADKSV VASSATSSLS YPSPVATSY PSPVTTSYPS

image



Background

Alternative Names	Early growth response protein 1 (EGR-1) (AT225) (Nerve growth factor-induced protein A) (NGFI-A) (Transcription factor ETR103) (Transcription factor Zif268) (Zinc finger protein 225) (Zinc finger protein Krox-24)
--------------------------	--

Note

For Research Use Only , Not for Diagnostic Use