


Product Name: Recombinant Human TAF1 Protein
Catalog #: PEH4804



Summary

Name	TAF1_HUMAN
Purity	Greater than 85% as determined by reducing SDS-PAGE
Endotoxin level	Please contact with the lab for this information
Tag	C-His
Accession #	P21675
Host	E.coli
Species	Human
Molecular Mass	19 kDa
Buffer	Lyophilized from Phosphate buffered saline
Shipping	The product is shipped on dry ice/polar packs. Upon receipt, store it immediately at the temperature listed below.
Stability&Storage	Store at $\leq -70^{\circ}\text{C}$, stable for 12 months after receipt. Store at $\leq -70^{\circ}\text{C}$, stable for 3 months under sterile conditions after opening. Please minimize freeze-thaw cycles.
Sequence	PQPVKKTVTGT TDADLRRLSL KNAKQLLRKF GVPEEEIKKL SRWEVIDVVR TMSTEQARSG EGPMSKFARG SRFVSAEHQE RYKEECQRIF DLQNKVLSST EVLSTDTDSS SAEDSDFEEM GKNIENMLQN KKTSSQLSRE REEQERKELQ

image



Background

Alternative Names	Transcription initiation factor TFIID subunit 1 (EC 2.3.1.48) (EC 2.7.11.1) (Cell cycle gene 1 protein) (TBP-associated factor 250 kDa) (p250) (Transcription initiation factor TFIID 250 kDa subunit) (TAF(II)250) (TAFII-250) (TAFII250)
--------------------------	--

Note

For Research Use Only , Not for Diagnostic Use

Product Name: Recombinant Human TAF1 Protein
Catalog #: PEH4804

