

Summary

| | |
|------------------------------|--|
| Name | ZBT25_HUMAN |
| Purity | Greater than 85% as determined by reducing SDS-PAGE |
| Endotoxin level | Please contact with the lab for this information |
| Tag | C-His |
| Accession # | P24278 |
| Host | E.coli |
| Species | Human |
| Molecular Mass | 23 kDa |
| Buffer | Phosphate buffered saline |
| Shipping | The product is shipped on dry ice/polar packs. Upon receipt, store it immediately at the temperature listed below. |
| Stability&Storage | Store at $\leq -70^{\circ}\text{C}$, stable for 6 months after receipt. Store at $\leq -70^{\circ}\text{C}$, stable for 3 months under sterile conditions after opening. Please minimize freeze-thaw cycles. |
| Sequence | VKQGLEVKEA PSSNSGNRAA VQGDHPQLQL SLAIGLDDGT ADQQRACPAT QALEEHQKPP VSIKQERCDP ESVISQSHPS PSSEVTGPTF TENSVKIHLC HYCGERFDSR SNLRQHLHTH VSGSLPFGVP ASILESNDLG EVHPLNENSE ALECRRLSSF IVKENEQQPD HTNRGTTEPL QISQVSLISK DTEPVELNCN |

image



Background

| | |
|--------------------------|--|
| Alternative Names | Zinc finger and BTB domain-containing protein 25 (Zinc finger protein 46) (Zinc finger protein KUP) |
|--------------------------|--|

Note

For Research Use Only , Not for Diagnostic Use

Product Name: Recombinant Human ZBTB25 Protein
Catalog #: PEH4808

