


Product Name: Recombinant Human ARID4A Protein
Catalog #: PEH4812



Summary

Name	ARI4A_HUMAN
Purity	Greater than 85% as determined by reducing SDS-PAGE
Endotoxin level	Please contact with the lab for this information
Tag	C-His
Accession #	P29374
Host	E.coli
Species	Human
Molecular Mass	56 kDa
Buffer	Phosphate buffered saline
Shipping	The product is shipped on dry ice/polar packs. Upon receipt, store it immediately at the temperature listed below.
Stability&Storage	Store at $\leq -70^{\circ}\text{C}$, stable for 6 months after receipt. Store at $\leq -70^{\circ}\text{C}$, stable for 3 months under sterile conditions after opening. Please minimize freeze-thaw cycles.
Sequence	TVGTDVSAKY RGAFCEAKIK TVKRLVKVKV LLKQDNTTQL VQDDQVKGPL RVGAIVETRT SDGSFQEAIL SKLTDASWYT VVFDGDGERT LRRTSLCLKG ERHFAESETL DQLPLTNPEH FGTPVIAKKT NRGRRSSLPV TEDEKEEESS EEDEDEKRRLL NDELLGKVVS VVSATERTEW YPALVISPC NDDITVKKDQ CLVRSFIDSK FYSIARKDIK EVDILNLPES ELSTKPGLQK ASIFLKTRVV PDNWKMDISE ILESSSSDDE DGPAEENDEE KEKEAKKTEE EVPEEELDPE ERDNFLQQLY KFMEDRGTPI NKPPVLGYKD LNLFLKFLRV YHGGGCDNID SGAVWKQIYM DLGIPILNSA ASYNVKTAYR KYLYGFEEYC RSANIQFRTV HHHEPKVKEE KKDLEESMEE ALKLDQEMPL TEVKSEPEEN IDNSSESERE EIELKSPRGR RRIARDVNSI

image



Background

Product Name: Recombinant Human ARID4A Protein
Catalog #: PEH4812



Alternative Names

AT-rich interactive domain-containing protein 4A (ARID domain-containing protein 4A) (Retinoblastoma-binding protein 1) (RBBP-1)

Note

For Research Use Only , Not for Diagnostic Use