


Product Name: Recombinant Human GTF2H1 Protein
Catalog #: PEH4819



Summary

Name	TF2H1_HUMAN
Purity	Greater than 85% as determined by reducing SDS-PAGE
Endotoxin level	Please contact with the lab for this information
Tag	C-His
Accession #	P32780
Host	E.coli
Species	Human
Molecular Mass	15 kDa
Buffer	Phosphate buffered saline
Shipping	The product is shipped on dry ice/polar packs. Upon receipt, store it immediately at the temperature listed below.
Stability&Storage	Store at $\leq -70^{\circ}\text{C}$, stable for 6 months after receipt. Store at $\leq -70^{\circ}\text{C}$, stable for 3 months under sterile conditions after opening. Please minimize freeze-thaw cycles.
Sequence	DLFAECAKID EKGLKTMVSL GVKNPPLDLT ALEDKPLDEG YGISSVPSAS NSKSIKENS N AAIKRNFNHH SAMVLAAGLR KQEAQNEQTS EPSNMDGNSG DADCFQPAVK RAKLQESIEY

image



Background

Alternative Names	General transcription factor IIH subunit 1 (Basic transcription factor 2 62 kDa subunit) (BTF2 p62) (General transcription factor IIH polypeptide 1) (TFIIH basal transcription factor complex p62 subunit)
--------------------------	---

Note

For Research Use Only , Not for Diagnostic Use

Product Name: Recombinant Human GTF2H1 Protein
Catalog #: PEH4819

