


Product Name: Recombinant Human RFC1 Protein
Catalog #: PEH4822



Summary

Name	RFC1_HUMAN
Purity	Greater than 85% as determined by reducing SDS-PAGE
Endotoxin level	Please contact with the lab for this information
Tag	C-His
Accession #	P35251
Host	E.coli
Species	Human
Molecular Mass	14 kDa
Buffer	Phosphate buffered saline
Shipping	The product is shipped on dry ice/polar packs. Upon receipt, store it immediately at the temperature listed below.
Stability&Storage	Store at $\leq -70^{\circ}\text{C}$, stable for 6 months after receipt. Store at $\leq -70^{\circ}\text{C}$, stable for 3 months under sterile conditions after opening. Please minimize freeze-thaw cycles.
Sequence	MVASKRKELS QNTDESLND EAIKQLQLD EDAELERQLH EDEEFARTLA MLDEEPTKK ARKDTEAGET FSSVQANLSK AEKHYPHKV KTAQVSDERK SYSPRKQSKY

image



Background

Alternative Names	Replication factor C subunit 1 (Activator 1 140 kDa subunit) (A1 140 kDa subunit) (Activator 1 large subunit) (Activator 1 subunit 1) (DNA-binding protein PO-GA) (Replication factor C 140 kDa subunit) (RF-C 140 kDa subunit) (RFC140) (Replication factor C large subunit)
--------------------------	---

Note

Product Name: Recombinant Human RFC1 Protein
Catalog #: PEH4822



For Research Use Only , Not for Diagnostic Use