


Product Name: Recombinant Human ZNF93 Protein
Catalog #: PEH4823



Summary

Name	ZNF93_HUMAN
Purity	Greater than 85% as determined by reducing SDS-PAGE
Endotoxin level	Please contact with the lab for this information
Tag	C-His
Accession #	P35789
Host	E.coli
Species	Human
Molecular Mass	16 kDa
Buffer	Phosphate buffered saline
Shipping	The product is shipped on dry ice/polar packs. Upon receipt, store it immediately at the temperature listed below.
Stability&Storage	Store at $\leq -70^{\circ}\text{C}$, stable for 6 months after receipt. Store at $\leq -70^{\circ}\text{C}$, stable for 3 months under sterile conditions after opening. Please minimize freeze-thaw cycles.
Sequence	ECGKAFKYSS TLSSHKRSHT GEKPYKCEEC GKAFVASSTL SKHEIHTGK KPYKCEECGK AFNQSSSLTK HKKIHTGEKP YKCEECGKAF NQSSSLTKHK KIHTGEKPYK CEECGKAFNQ

image



Background

Alternative Names Zinc finger protein 93 (Zinc finger protein 505) (Zinc finger protein HTF34)

Note

For Research Use Only , Not for Diagnostic Use