


Product Name: Recombinant Human ZNF76 Protein
Catalog #: PEH4824



Summary

Name	ZNF76_HUMAN
Purity	Greater than 85% as determined by reducing SDS-PAGE
Endotoxin level	Please contact with the lab for this information
Tag	C-His
Accession #	P36508
Host	E.coli
Species	Human
Molecular Mass	16 kDa
Buffer	Phosphate buffered saline
Shipping	The product is shipped on dry ice/polar packs. Upon receipt, store it immediately at the temperature listed below.
Stability&Storage	Store at $\leq -70^{\circ}\text{C}$, stable for 6 months after receipt. Store at $\leq -70^{\circ}\text{C}$, stable for 3 months under sterile conditions after opening. Please minimize freeze-thaw cycles.
Sequence	GQVIQLEDGT TAYIHQVTVQ KEALSFEDGQ PVQLEDGSMA YIHRTPREGY DPSTLEAVQL EDGSTAYIHH PVAVPSESTI LAVQTEVGLD DLAAEDDEGF SADAVVALEQ YASKVLHDSQ IPRNGKGQQV

image



Background

Alternative Names Zinc finger protein 76 (Zinc finger protein 523)

Note

For Research Use Only , Not for Diagnostic Use