


Product Name: Recombinant Human ADCY8 Protein
Catalog #: PEH4825



Summary

Name	ADCY8_HUMAN
Purity	Greater than 85% as determined by reducing SDS-PAGE
Endotoxin level	Please contact with the lab for this information
Tag	C-His
Accession #	P40145
Host	E.coli
Species	Human
Molecular Mass	18 kDa
Buffer	Lyophilized from Phosphate buffered saline
Shipping	The product is shipped on dry ice/polar packs. Upon receipt, store it immediately at the temperature listed below.
Stability&Storage	Store at $\leq -70^{\circ}\text{C}$, stable for 12 months after receipt. Store at $\leq -70^{\circ}\text{C}$, stable for 3 months under sterile conditions after opening. Please minimize freeze-thaw cycles.
Sequence	KELREHNENM LRNILPSHVA RHFLEKDRDN EELYSQSYDA VGVMFASIPG FADFYSQTEM NNQGVECLRL LNEIIADFDE LLGEDRFQDI EKIKTIGSTY MAVSGLSPEK QQCEDKKGHL CALADFSLAL

image



Background

Alternative Names	Adenylate cyclase type 8 (EC 4.6.1.1) (ATP pyrophosphate-lyase 8) (Adenylate cyclase type VIII) (Adenylyl cyclase 8) (AC8) (Ca ²⁺)/calmodulin-activated adenylyl cyclase)
--------------------------	---

Note

For Research Use Only , Not for Diagnostic Use

Product Name: Recombinant Human ADCY8 Protein
Catalog #: PEH4825

