


Product Name: Recombinant Human PBX3 Protein
Catalog #: PEH4826



Summary

Name	PBX3_HUMAN
Purity	Greater than 85% as determined by reducing SDS-PAGE
Endotoxin level	Please contact with the lab for this information
Tag	C-His
Accession #	P40426
Host	E.coli
Species	Human
Molecular Mass	19 kDa
Buffer	Lyophilized from Phosphate buffered saline
Shipping	The product is shipped on dry ice/polar packs. Upon receipt, store it immediately at the temperature listed below.
Stability&Storage	Store at $\leq -70^{\circ}\text{C}$, stable for 12 months after receipt. Store at $\leq -70^{\circ}\text{C}$, stable for 3 months under sterile conditions after opening. Please minimize freeze-thaw cycles.
Sequence	NFSKQATEIL NEYFYSHLSN PYPSEEAKEE LAKKCSITVS QVSNWFGNKR IRYKKNIGKF QEEANLYAAK TAVTAAHAVA AAVQNNQTNS PTPNSGSSG SFNLPNSGDM FMNMQSLNGD SYQGSQVGAN VQSQVDTLRH VINQTGGYSD

image



Background

Alternative Names Pre-B-cell leukemia transcription factor 3 (Homeobox protein PBX3)

Note

For Research Use Only , Not for Diagnostic Use