


Product Name: Recombinant Human HELZ Protein
Catalog #: PEH4827



Summary

Name	HELZ_HUMAN
Purity	Greater than 85% as determined by reducing SDS-PAGE
Endotoxin level	Please contact with the lab for this information
Tag	C-His
Accession #	P42694
Host	E.coli
Species	Human
Molecular Mass	14 kDa
Buffer	Phosphate buffered saline
Shipping	The product is shipped on dry ice/polar packs. Upon receipt, store it immediately at the temperature listed below.
Stability&Storage	Store at $\leq -70^{\circ}\text{C}$, stable for 6 months after receipt. Store at $\leq -70^{\circ}\text{C}$, stable for 3 months under sterile conditions after opening. Please minimize freeze-thaw cycles.
Sequence	AVDPRIITHQ AAMAYNMNLL QTHGRGSPIP YGLGHHPPVT IGQPQNQHQE KDQHEQNRNG KSDTNNSGPE INKIRTPEKK PTEPKQVDLE SNPQNRSPES RPSVVYPSTK

image



Background

Alternative Names Probable helicase with zinc finger domain (EC 3.6.4.-) (Down-regulated in human cancers protein)

Note

For Research Use Only , Not for Diagnostic Use