


Product Name: Recombinant Human RBM34 Protein
Catalog #: PEH4828



Summary

Name	RBM34_HUMAN
Purity	Greater than 85% as determined by reducing SDS-PAGE
Endotoxin level	Please contact with the lab for this information
Tag	C-His
Accession #	P42696
Host	E.coli
Species	Human
Molecular Mass	25 kDa
Buffer	Lyophilized from Phosphate buffered saline
Shipping	The product is shipped on dry ice/polar packs.Upon receipt, store it immediately at the temperature listed below.
Stability&Storage	Store at $\leq -70^{\circ}\text{C}$, stable for 12 months after receipt.Store at $\leq -70^{\circ}\text{C}$, stable for 3 months under sterile conditions after opening. Please minimize freeze-thaw cycles.
Sequence	KTKRNEEEES TSQIERPLSQ EPAKKVKAKK KHTNAEKKLA DRESALASAD LEEEIHQKQG QKRKNSQPGV KVADRKILDD TEDTVVSQRK KIQINQEEER LKNERTVFIG NLPVTCNKKK LKSFFKEYGQ IESVRFRLSI PAEGTLSKKL AAIKRKIHPD QKNINAYVVF KEESAATQAL KRNGAQIADG FRIRVDLASE

image



Background

Alternative Names RNA-binding protein 34 (RNA-binding motif protein 34)

Note

For Research Use Only , Not for Diagnostic Use

Product Name: Recombinant Human RBM34 Protein
Catalog #: PEH4828

