


**Product Name: Recombinant Human GUCA1A Protein**  
**Catalog #: PEH4829**



## Summary

|                              |   |
|------------------------------|---|
| <b>Name</b>                  | GUC1A_HUMAN   |
| <b>Purity</b>                | Greater than 85% as determined by reducing SDS-PAGE   |
| <b>Endotoxin level</b>       | Please contact with the lab for this information  |
| <b>Tag</b>                   | C-His   |
| <b>Accession #</b>           | P43080  |
| <b>Host</b>                  | E.coli  |
| <b>Species</b>               | Human   |
| <b>Molecular Mass</b>        | 24 kDa  |
| <b>Buffer</b>                | Phosphate buffered saline   |
| <b>Shipping</b>              | The product is shipped on dry ice/polar packs. Upon receipt, store it immediately at the temperature listed below.  |
| <b>Stability&amp;Storage</b> | Store at $\leq -70^{\circ}\text{C}$ , stable for 6 months after receipt. Store at $\leq -70^{\circ}\text{C}$ , stable for 3 months under sterile conditions after opening. Please minimize freeze-thaw cycles.                    |
| <b>Sequence</b>              | MGNVMEGKSV EELSSTECHQ WYKKFMTECP SGQLTYEFR QFFGLKNLSP<br>SASQYVEQMF ETDFNKGDY IDMEYVAAL SLVLKGKVEQ KLRWYFKLYD<br>VDGNGCIDRD ELLTIQAIR AINPCSDTTM TAEFTDTVF SKIDVNGDGE<br>LSLEEFIEGV QKDQMLLDTL TRSLDLTRIV RRLQNGEQDE EGADEAAEEA G |

## image



## Background

**Alternative Names** Guanylyl cyclase-activating protein 1 (GCAP 1) (Guanylate cyclase activator 1A)

## Note

For Research Use Only , Not for Diagnostic Use

**Product Name: Recombinant Human GUCA1A Protein**  
**Catalog #: PEH4829**

---

