


Product Name: Recombinant Human LHX1 Protein
Catalog #: PEH4832



Summary

Name	LHX1_HUMAN
Purity	Greater than 85% as determined by reducing SDS-PAGE
Endotoxin level	Please contact with the lab for this information
Tag	C-His
Accession #	P48742
Host	E.coli
Species	Human
Molecular Mass	18 kDa
Buffer	Phosphate buffered saline
Shipping	The product is shipped on dry ice/polar packs. Upon receipt, store it immediately at the temperature listed below.
Stability&Storage	Store at $\leq -70^{\circ}\text{C}$, stable for 6 months after receipt. Store at $\leq -70^{\circ}\text{C}$, stable for 3 months under sterile conditions after opening. Please minimize freeze-thaw cycles.
Sequence	HSATTGSDPS LSPDSQDPSQ DDAKDESAN VSDKEAGSNE NDDQNLGAKR RGPRTTIKAK QLETLKAAFA ATPKPTRHIR EQLAQETGLN MRVIQVWFQN RRSKERRMKQ LALGARRHA FFRSPRRMRP LVDRLEPGEL

image



Background

Alternative Names LIM/homeobox protein Lhx1 (LIM homeobox protein 1) (Homeobox protein Lim-1) (hLim-1)

Note

For Research Use Only , Not for Diagnostic Use