


Product Name: Recombinant Human TUB Protein
Catalog #: PEH4836



Summary

Name	TUB_HUMAN
Purity	Greater than 85% as determined by reducing SDS-PAGE
Endotoxin level	Please contact with the lab for this information
Tag	C-His
Accession #	P50607
Host	E.coli
Species	Human
Molecular Mass	22 kDa
Buffer	Phosphate buffered saline
Shipping	The product is shipped on dry ice/polar packs. Upon receipt, store it immediately at the temperature listed below.
Stability&Storage	Store at $\leq -70^{\circ}\text{C}$, stable for 6 months after receipt. Store at $\leq -70^{\circ}\text{C}$, stable for 3 months under sterile conditions after opening. Please minimize freeze-thaw cycles.
Sequence	NLRQQKLDLRQ RALLEQKQKK KRQEPLMVQA NADGRPRSRR ARQSEEQAPL VESYLSSSGS TSYQVQEADS LASVQLGATR PTAPASAKRT KAAATAGGQG GAARKEKKGK HKGTSGPAAL AEDKSEAQGP VQILTVGQSD HAQDAGETAA GGGERPSGQD LRATMQRKGI SSSMSFDEDE EDEEENSSSS

image



Background

Alternative Names Tubby protein homolog

Note

For Research Use Only , Not for Diagnostic Use

Product Name: Recombinant Human TUB Protein
Catalog #: PEH4836

