


Product Name: Recombinant Human GBX2 Protein
Catalog #: PEH4841



Summary

Name	GBX2_HUMAN
Purity	Greater than 85% as determined by reducing SDS-PAGE
Endotoxin level	Please contact with the lab for this information
Tag	C-His
Accession #	P52951
Host	E.coli
Species	Human
Molecular Mass	18 kDa
Buffer	Phosphate buffered saline
Shipping	The product is shipped on dry ice/polar packs. Upon receipt, store it immediately at the temperature listed below.
Stability&Storage	Store at $\leq -70^{\circ}\text{C}$, stable for 6 months after receipt. Store at $\leq -70^{\circ}\text{C}$, stable for 3 months under sterile conditions after opening. Please minimize freeze-thaw cycles.
Sequence	LESVDYSSD DNLTGQAAHK EEDPGHALEE TPPSSGAAGS TTSTGKNRRR RTAFTSEQLL ELEKEFHCKK YLSLTERSQI AHALKLSEVQ VKIWFQNRRA KWKRVKAGNA NSKTGEP SRN PKIVVPIPVH VSRFAIRSQH

image



Background

Alternative Names Homeobox protein GBX-2 (Gastrulation and brain-specific homeobox protein 2)

Note

For Research Use Only , Not for Diagnostic Use