

**Product Name: Recombinant Human PTPN5 Protein**  
**Catalog #: PEH4846**



## Summary

<b>Name</b>	PTN5_HUMAN
<b>Purity</b>	Greater than 85% as determined by reducing SDS-PAGE
<b>Endotoxin level</b>	Please contact with the lab for this information
<b>Tag</b>	C-His
<b>Accession #</b>	P54829
<b>Host</b>	E.coli
<b>Species</b>	Human
<b>Molecular Mass</b>	15 kDa
<b>Buffer</b>	Lyophilized from Phosphate buffered saline
<b>Shipping</b>	The product is shipped on dry ice/polar packs. Upon receipt, store it immediately at the temperature listed below.
<b>Stability&amp;Storage</b>	Store at $\leq -70^{\circ}\text{C}$ , stable for 12 months after receipt. Store at $\leq -70^{\circ}\text{C}$ , stable for 3 months under sterile conditions after opening. Please minimize freeze-thaw cycles.
<b>Sequence</b>	RRQSVSRQPS FTYSEWMEEK IEDDFLDLDP VPETPVFDCV MDIKPEADPT SLTVKSMGLQ ERRGSNVSLT LDMCTPGCNE EFGYLMSPR EESAREYLLS ASRVLQAEEL

## image



## Background

<b>Alternative Names</b>	Tyrosine-protein phosphatase non-receptor type 5 (EC 3.1.3.48) (Neural-specific protein-tyrosine phosphatase) (Striatum-enriched protein-tyrosine phosphatase) (STEP)
--------------------------	---

## Note

For Research Use Only , Not for Diagnostic Use

**Product Name: Recombinant Human PTPN5 Protein**  
**Catalog #: PEH4846**

---

