


Product Name: Recombinant Human PRDM15 Protein
Catalog #: PEH4849



Summary

Name	PRD15_HUMAN
Purity	Greater than 85% as determined by reducing SDS-PAGE
Endotoxin level	Please contact with the lab for this information
Tag	C-His
Accession #	P57071
Host	E.coli
Species	Human
Molecular Mass	20 kDa
Buffer	Phosphate buffered saline
Shipping	The product is shipped on dry ice/polar packs. Upon receipt, store it immediately at the temperature listed below.
Stability&Storage	Store at $\leq -70^{\circ}\text{C}$, stable for 6 months after receipt. Store at $\leq -70^{\circ}\text{C}$, stable for 3 months under sterile conditions after opening. Please minimize freeze-thaw cycles.
Sequence	VLAVRIDDL LHPETTTIDA SSIGIVQPEL TLEQEDLAEG KHGKAAKRSH KRKQKPEEEA GAPVPEDATF SEYSEKETEF TGSVGDETNS AVQSIQQVVV TLGDPNVTTP SSSVGLTNIT VTPITTAAT QFTNLQPVAV GHLTTPERQL QLDNSILTVT FDTVSGSAML

image



Background

Alternative Names	PR domain zinc finger protein 15 (EC 2.1.1.-) (PR domain-containing protein 15) (Zinc finger protein 298)
--------------------------	---

Note

For Research Use Only , Not for Diagnostic Use

Product Name: Recombinant Human PRDM15 Protein
Catalog #: PEH4849

