


Product Name: Recombinant Human EVC Protein
Catalog #: PEH4850



Summary

Name	EVC_HUMAN
Purity	Greater than 85% as determined by reducing SDS-PAGE
Endotoxin level	Please contact with the lab for this information
Tag	C-His
Accession #	P57679
Host	E.coli
Species	Human
Molecular Mass	37 kDa
Buffer	Phosphate buffered saline
Shipping	The product is shipped on dry ice/polar packs. Upon receipt, store it immediately at the temperature listed below.
Stability&Storage	Store at $\leq -70^{\circ}\text{C}$, stable for 6 months after receipt. Store at $\leq -70^{\circ}\text{C}$, stable for 3 months under sterile conditions after opening. Please minimize freeze-thaw cycles.
Sequence	KKSDDELYQK ILSKQEKDLE ELEKGLQVKL SNTEMSGAGD SEYITLADVE KKEREYSEQL IDNMEAFWKQ MANIQHFLVD QFKCSSSKAR QLMMTLTERM IAAEGLLCDS QELQALDALE RTMGRAHMAK VIEFLKLQVQ EETRCRLAAI SHGLELLAGE GKLSGRQKEE LLTQQHKAFW QEAERFSREF VQRGKDLVTA SLAHQVEGTA KLTLAQEEEQ RSFLAEAQPT ADPEKFLEAF HEVLERQRLM QCDLEEEENV RATEAVVALC QELYFSTVDT FQKFVDALFL QTLPGMTGLP

image



Background

Alternative Names Ellis-van Creveld syndrome protein (DWF-1)

Product Name: Recombinant Human EVC Protein
Catalog #: PEH4850



Note

For Research Use Only , Not for Diagnostic Use