

Product Name: Recombinant Human GJA9 Protein
Catalog #: PEH4852



Summary

Name	CXA9_HUMAN
Purity	Greater than 85% as determined by reducing SDS-PAGE
Endotoxin level	Please contact with the lab for this information
Tag	C-His
Accession #	P57773
Host	E.coli
Species	Human
Molecular Mass	21 kDa
Buffer	Phosphate buffered saline
Shipping	The product is shipped on dry ice/polar packs. Upon receipt, store it immediately at the temperature listed below.
Stability&Storage	Store at $\leq -70^{\circ}\text{C}$, stable for 6 months after receipt. Store at $\leq -70^{\circ}\text{C}$, stable for 3 months under sterile conditions after opening. Please minimize freeze-thaw cycles.
Sequence	GLWGKYKLLK EHNEFHANKA KQNVAKYQST SANSLKRLPS APDYNLLVEK QTHTAVYPSL NSSSVFQPNP DNHSVNDEKC ILDEQETVLS NEISTLSTSC SHFQHISSNN NKDTHKIFGK ELNGNQLMEK RETEGKDSKR NYYSRGRHSI PGVAIDGENN MRQSPQTVFS

image



Background

Alternative Names Gap junction alpha-9 protein (Connexin-58) (Cx58) (Connexin-59) (Cx59) (Gap junction alpha-10 protein)

Note

For Research Use Only , Not for Diagnostic Use

Product Name: Recombinant Human GJA9 Protein
Catalog #: PEH4852

