


Product Name: Recombinant Human ZNF280A Protein
Catalog #: PEH4853



Summary

Name	Z280A_HUMAN
Purity	Greater than 85% as determined by reducing SDS-PAGE
Endotoxin level	Please contact with the lab for this information
Tag	C-His
Accession #	P59817
Host	E.coli
Species	Human
Molecular Mass	27 kDa
Buffer	Phosphate buffered saline
Shipping	The product is shipped on dry ice/polar packs. Upon receipt, store it immediately at the temperature listed below.
Stability&Storage	Store at $\leq -70^{\circ}\text{C}$, stable for 6 months after receipt. Store at $\leq -70^{\circ}\text{C}$, stable for 3 months under sterile conditions after opening. Please minimize freeze-thaw cycles.
Sequence	ESPKKNLRES KQREEDDEDP DLIYVGVEHV HRDAEVLVFG MISNSKPVVS NILNRVTPGS NSRRKKGHFR QYPAHVSQPA NHVTSMAKAI MPVSLSEGRS TDSPVTMKSS SEPGYKMSSP QVVSPSSSDS LPPGTQCLVG AMVSGGGRNE SSPDSKRLST SDINSRDSKR VKLRDGIPGV PSLAVVPSDM SSTISTNTPS QGICNSSNHV QNGVTFPWPD ANGKAHFNL

image



Background

Alternative Names	Zinc finger protein 280A (3'OY11.1) (Suppressor of hairy wing homolog 1) (Zinc finger protein 636)
--------------------------	--

Note

Product Name: Recombinant Human ZNF280A Protein
Catalog #: PEH4853



For Research Use Only , Not for Diagnostic Use