


Product Name: Recombinant Human ERVK-8 Protein
Catalog #: PEH4856



Summary

Name	GAK8_HUMAN
Purity	Greater than 85% as determined by reducing SDS-PAGE
Endotoxin level	Please contact with the lab for this information
Tag	C-His
Accession #	P62685
Host	E.coli
Species	Human
Molecular Mass	21 kDa
Buffer	Phosphate buffered saline
Shipping	The product is shipped on dry ice/polar packs. Upon receipt, store it immediately at the temperature listed below.
Stability&Storage	Store at $\leq -70^{\circ}\text{C}$, stable for 6 months after receipt. Store at $\leq -70^{\circ}\text{C}$, stable for 3 months under sterile conditions after opening. Please minimize freeze-thaw cycles.
Sequence	KKSQKETESL HCEYVAEPVM AQSTQNV DYN QLQEV IYPET LKLEGKGPEL VGPSKPRG TSPLAGQVP VTLQPQKQVK ENKTQPPVAY QYWPPAELQY RPPESQYGY PGMPPAPQGR EPYPQPPTRR LNPTAPPSRQ GSELHEIDK SRKEGDTEAW QFPVTLEPMP PGEQAQEGEP

image



Background

Alternative Names	Endogenous retrovirus group K member 8 Gag polyprotein (HERV-K115 Gag protein) (HERV-K_8p23.1 provirus ancestral Gag polyprotein) (Gag polyprotein)
--------------------------	---

Note

Product Name: Recombinant Human ERVK-8 Protein
Catalog #: PEH4856



For Research Use Only , Not for Diagnostic Use