


Product Name: Recombinant Human ERVK-9 Protein
Catalog #: PEH4859



Summary

Name	GAK9_HUMAN
Purity	Greater than 85% as determined by reducing SDS-PAGE
Endotoxin level	Please contact with the lab for this information
Tag	C-His
Accession #	P63126
Host	E.coli
Species	Human
Molecular Mass	23 kDa
Buffer	Phosphate buffered saline
Shipping	The product is shipped on dry ice/polar packs. Upon receipt, store it immediately at the temperature listed below.
Stability&Storage	Store at $\leq -70^{\circ}\text{C}$, stable for 6 months after receipt. Store at $\leq -70^{\circ}\text{C}$, stable for 3 months under sterile conditions after opening. Please minimize freeze-thaw cycles.
Sequence	KKSQKETESL HCEYVAEPVM AQSTQNV DYN QLQEV IYPET LKLEGKGPEL VGPSKPRG TSPLAGQVP VTLQPQKQVK ENKTQPPVAY QYWPPAELQY RPPESQYGY PGMPPAPQGR APYPQPPTRR LNPTAPPSRQ GSELHEIDK SRKEGDTEAW QFPVTLEPMP PGEQAQEGEP PTVEARYKSF

image



Background

Alternative Names	Endogenous retrovirus group K member 9 Gag polyprotein (HERV-K(C6) Gag protein) (HERV-K109 Gag protein) (HERV-K_6q14.1 provirus ancestral Gag polyprotein) (Gag polyprotein)
--------------------------	--

Note

Product Name: Recombinant Human ERVK-9 Protein
Catalog #: PEH4859



For Research Use Only , Not for Diagnostic Use