


Product Name: Recombinant Human ERVK-7 Protein
Catalog #: PEH4861



Summary

Name	GAK7_HUMAN
Purity	Greater than 85% as determined by reducing SDS-PAGE
Endotoxin level	Please contact with the lab for this information
Tag	C-His
Accession #	P63130
Host	E.coli
Species	Human
Molecular Mass	21 kDa
Buffer	Lyophilized from Phosphate buffered saline
Shipping	The product is shipped on dry ice/polar packs.Upon receipt, store it immediately at the temperature listed below.
Stability&Storage	Store at $\leq -70^{\circ}\text{C}$, stable for 12 months after receipt.Store at $\leq -70^{\circ}\text{C}$, stable for 3 months under sterile conditions after opening. Please minimize freeze-thaw cycles.
Sequence	KKSQKETEG L HCEYVAEPVM AQSTQNV DYN QLQEV IYPET LKLEGKGPEL VGPSSESKPRG TSHLPAGQVP VTLQPQKQVK ENKTQPPVAY QYWPPAELQY RPPPE SQYGY PGMPPAPQGR APYPQP PTRR LNPTAPPSRQ GSELHEIDK SRKEGDTEAW QFPVTLEPMP

image



Background

Alternative Names	Endogenous retrovirus group K member 7 Gag polyprotein (HERV-K(III) Gag protein) (HERV-K102 Gag protein) (HERV-K_1q22 provirus ancestral Gag polyprotein) (Gag polyprotein)
--------------------------	---

Note

Product Name: Recombinant Human ERVK-7 Protein
Catalog #: PEH4861



For Research Use Only , Not for Diagnostic Use