

Product Name: Recombinant Human KRT83 Protein
Catalog #: PEH4867



Summary

Name	KRT83_HUMAN
Purity	Greater than 85% as determined by reducing SDS-PAGE
Endotoxin level	Please contact with the lab for this information
Tag	C-His
Accession #	P78385
Host	E.coli
Species	Human
Molecular Mass	40 kDa
Buffer	Phosphate buffered saline
Shipping	The product is shipped on dry ice/polar packs. Upon receipt, store it immediately at the temperature listed below.
Stability&Storage	Store at $\leq -70^{\circ}\text{C}$, stable for 6 months after receipt. Store at $\leq -70^{\circ}\text{C}$, stable for 3 months under sterile conditions after opening. Please minimize freeze-thaw cycles.
Sequence	SLLTPLNLEI DPNAQCVKQE EKEQIKSLNS RFAAFIDKVR FLEQQNKLL TKLQFYQNRE CCQSNLEPLF AGYIETLRRE AECVEADSGR LASELNHVQE VLEGYKKKYE EEVALRATAE NEFVALKKDV DCAYLRKSDL EANVEALIQE IDFLRRLYEE EIRILQSHIS DTSVVVKLDN SRDLNMDCIV AEIKAQYDDI ATRSRAEAE WYRSKCEEMK ATVIRHGETL RRTKEEINEL NRMIQRLTAE VENAKCQNSK LEAAVAQSEQ QGEAALSDAR CKLAELEGAL QKAKQDMACL IREYQEVMSN KLGLDIEIAT YRRLLEGEEQ

image



Background

Alternative Names	Keratin, type II cuticular Hb3 (Hair keratin K2.10) (Keratin-83) (K83) (Type II hair keratin Hb3) (Type-II keratin Kb23)
--------------------------	--

Product Name: Recombinant Human KRT83 Protein
Catalog #: PEH4867



Note

For Research Use Only , Not for Diagnostic Use