


Product Name: Recombinant Human ADARB1 Protein
Catalog #: PEH4869



Summary

Name	RED1_HUMAN
Purity	Greater than 85% as determined by reducing SDS-PAGE
Endotoxin level	Please contact with the lab for this information
Tag	C-His
Accession #	P78563
Host	E.coli
Species	Human
Molecular Mass	22 kDa
Buffer	Phosphate buffered saline
Shipping	The product is shipped on dry ice/polar packs. Upon receipt, store it immediately at the temperature listed below.
Stability&Storage	Store at $\leq -70^{\circ}\text{C}$, stable for 6 months after receipt. Store at $\leq -70^{\circ}\text{C}$, stable for 3 months under sterile conditions after opening. Please minimize freeze-thaw cycles.
Sequence	SSSTDVKENR NLDNVSPKDG STPGPGEQSQ LSNGGGGGPG RKRPLEEGSN GHSKYRLKRR RKTGPVLPK NALMQLNEIK PGLQYTLLSQ TGPVHAPLFV MSVEVNGQVF EGS GPTKKKA KLHAAEKALR SFVQFPNASE AHLAMGRTLS VNTDFTSDQA DFPDTLFGF ETPDKAEPFF

image



Background

Alternative Names Double-stranded RNA-specific editase 1 (EC 3.5.4.37) (RNA-editing deaminase 1) (RNA-editing enzyme 1) (dsRNA adenosine deaminase)

Note

For Research Use Only , Not for Diagnostic Use

Product Name: Recombinant Human ADARB1 Protein
Catalog #: PEH4869

