


Product Name: Recombinant Human RBM10 Protein
Catalog #: PEH4873



Summary

Name	RBM10_HUMAN
Purity	Greater than 85% as determined by reducing SDS-PAGE
Endotoxin level	Please contact with the lab for this information
Tag	C-His
Accession #	P98175
Host	E.coli
Species	Human
Molecular Mass	20 kDa
Buffer	Phosphate buffered saline
Shipping	The product is shipped on dry ice/polar packs. Upon receipt, store it immediately at the temperature listed below.
Stability&Storage	Store at $\leq -70^{\circ}\text{C}$, stable for 6 months after receipt. Store at $\leq -70^{\circ}\text{C}$, stable for 3 months under sterile conditions after opening. Please minimize freeze-thaw cycles.
Sequence	QGESEEEEEED EEEEEKASNI VMLRMLPQAA TEDDIRGQLQ SHGVQAREVR LMRNKSSGQS RGFAFVEFSH LQDATRWMEA NQHSLNILGQ KVSMHYSDPK PKINEDWLCN KCGVQNFKRR EKCFKCGVPK SEAEQKLPLG TRLDQQTLPL

image



Background

Alternative Names RNA-binding protein 10 (G patch domain-containing protein 9) (RNA-binding motif protein 10) (RNA-binding protein S1-1) (S1-1)

Note

For Research Use Only , Not for Diagnostic Use