


Product Name: Recombinant Human ATP8B2 Protein
Catalog #: PEH4874



Summary

| | |
|------------------------------|---|
| Name | AT8B2_HUMAN |
| Purity | Greater than 85% as determined by reducing SDS-PAGE |
| Endotoxin level | Please contact with the lab for this information |
| Tag | C-His |
| Accession # | P98198 |
| Host | E.coli |
| Species | Human |
| Molecular Mass | 21 kDa |
| Buffer | Lyophilized from Phosphate buffered saline |
| Shipping | The product is shipped on dry ice/polar packs. Upon receipt, store it immediately at the temperature listed below. |
| Stability&Storage | Store at $\leq -70^{\circ}\text{C}$, stable for 12 months after receipt. Store at $\leq -70^{\circ}\text{C}$, stable for 3 months under sterile conditions after opening. Please minimize freeze-thaw cycles. |
| Sequence | DLDEEYEEW AERRLQASLA QDSREDRLAS IYEEVENMMM LLGATAIEDK LQQGVPETIA LLTLANIKIW VLTGDKQETA VNIGYSCKML TDDMTEVFIV TGHTVLEVRE ELRKAREKMM DSSRSVGNGF TYQDKLSSSK LTSVLEAVAG EYALVINGHS LAHALEADME |

image



Background

Alternative Names Phospholipid-transporting ATPase ID (EC 7.6.2.1) (ATPase class I type 8B member 2) (P4-ATPase flippase complex alpha subunit ATP8B2)

Note

For Research Use Only , Not for Diagnostic Use

Product Name: Recombinant Human ATP8B2 Protein
Catalog #: PEH4874

