



Product Name: Recombinant Human TLE3 Protein
Catalog #: PEH4883



Summary

Name	TLE3_HUMAN
Purity	Greater than 85% as determined by reducing SDS-PAGE
Endotoxin level	Please contact with the lab for this information
Tag	C-His
Accession #	Q04726
Host	E.coli
Species	Human
Molecular Mass	23 kDa
Buffer	Phosphate buffered saline
Shipping	The product is shipped on dry ice/polar packs.Upon receipt, store it immediately at the temperature listed below.
Stability&Storage	Store at $\leq -70^{\circ}\text{C}$, stable for 6 months after receipt.Store at $\leq -70^{\circ}\text{C}$, stable for 3 months under sterile conditions after opening. Please minimize freeze-thaw cycles.
Sequence	IAKRLNTILA QIMPFLSQEH QQQVAQAVER AKQVTMTELN AIIGQQQLQA QHLSHATHGP PVQLPPHPSPG LQPPGIPPVT GSSSGLLAG ALGSAHLTV KDEKNHHELD HRERESSANN SVSPSESLRA SEKHRSADY SMEAKKRKAE EKDSLRYDS DGDKSDDLTV DVSNEPATP RVSPAHSPPPE

image



Background

Alternative Names	Transducin-like enhancer protein 3 (Enhancer of split groucho-like protein 3) (ESG3)
--------------------------	--

Note

For Research Use Only , Not for Diagnostic Use

Product Name: Recombinant Human TLE3 Protein
Catalog #: PEH4883

