


Summary

Name	MTG8_HUMAN
Purity	Greater than 85% as determined by reducing SDS-PAGE
Endotoxin level	Please contact with the lab for this information
Tag	C-His
Accession #	Q06455
Host	E.coli
Species	Human
Molecular Mass	21 kDa
Buffer	Phosphate buffered saline
Shipping	The product is shipped on dry ice/polar packs. Upon receipt, store it immediately at the temperature listed below.
Stability&Storage	Store at $\leq -70^{\circ}\text{C}$, stable for 6 months after receipt. Store at $\leq -70^{\circ}\text{C}$, stable for 3 months under sterile conditions after opening. Please minimize freeze-thaw cycles.
Sequence	DLKKGSSSS SHSRQQSPVN PDPVALDAHR EFLHRPASGY VPEEIWKKA EAVNEVKRQA MTELQKAVSE AERKAHDMIT TERA KMERTV AEAKRQAAED ALAVINQQED SSESCWNCGR KASETCGCGN TARYCGSFCQ HKDWEKHHHI CGQTLQAQQQ GDTPAVSSSV TPNSGAGSPM

image



Background

Alternative Names	Protein CBFA2T1 (Cyclin-D-related protein) (Eighty one protein) (Protein ETO) (Protein MTG8) (Zinc finger MYND domain-containing protein 2)
--------------------------	---

Note

Product Name: Recombinant Human RUNX1T1 Protein
Catalog #: PEH4885



For Research Use Only , Not for Diagnostic Use