


Product Name: Recombinant Human ZNF33A Protein
Catalog #: PEH4886



Summary

Name	ZN33A_HUMAN
Purity	Greater than 85% as determined by reducing SDS-PAGE
Endotoxin level	Please contact with the lab for this information
Tag	C-His
Accession #	Q06730
Host	E.coli
Species	Human
Molecular Mass	16 kDa
Buffer	Phosphate buffered saline
Shipping	The product is shipped on dry ice/polar packs. Upon receipt, store it immediately at the temperature listed below.
Stability&Storage	Store at $\leq -70^{\circ}\text{C}$, stable for 6 months after receipt. Store at $\leq -70^{\circ}\text{C}$, stable for 3 months under sterile conditions after opening. Please minimize freeze-thaw cycles.
Sequence	HDETHTQEK N EVLKNRNTLS HHEETLQHEK IQTLEHNF EY SICQETLLEK AVFNTQKRE N AEENNCDYNE FGRTLCDSSS LLFHQISPSR DNHYEFS DCE KFLCVKSTLS KPHGVSMKHY

image



Background

Alternative Names	Zinc finger protein 33A (Zinc finger and ZAK-associated protein with KRAB domain) (ZZaPK) (Zinc finger protein 11A) (Zinc finger protein KOX31)
--------------------------	---

Note

For Research Use Only , Not for Diagnostic Use