


Product Name: Recombinant Human USP17L24 Protein
Catalog #: PEH4889



Summary

Name	U17LO_HUMAN
Purity	Greater than 85% as determined by reducing SDS-PAGE
Endotoxin level	Please contact with the lab for this information
Tag	C-His
Accession #	Q0WX57
Host	E.coli
Species	Human
Molecular Mass	17 kDa
Buffer	Lyophilized from Phosphate buffered saline
Shipping	The product is shipped on dry ice/polar packs. Upon receipt, store it immediately at the temperature listed below.
Stability&Storage	Store at $\leq -70^{\circ}\text{C}$, stable for 12 months after receipt. Store at $\leq -70^{\circ}\text{C}$, stable for 3 months under sterile conditions after opening. Please minimize freeze-thaw cycles.
Sequence	RRATQGELKR DHPCLQAPEL DEHLVERATQ ESTLDHWKFL QEQNKTKEF NVRKVEGTLP PDVLVIHQSK YKCGMKNHHP EQQSSLLNLS SSTPTHQESM NTGTLASLRG RARRSKGKNK

image



Background

Alternative Names	Ubiquitin carboxyl-terminal hydrolase 17-like protein 24 (EC 3.4.19.12) (Deubiquitinating enzyme 17) (Ubiquitin thioesterase 17) (Ubiquitin-specific-processing protease 17)
--------------------------	--

Note

For Research Use Only , Not for Diagnostic Use

Product Name: Recombinant Human USP17L24 Protein
Catalog #: PEH4889

