


## Summary

<b>Name</b>	SF3A3_HUMAN
<b>Purity</b>	Greater than 85% as determined by reducing SDS-PAGE
<b>Endotoxin level</b>	Please contact with the lab for this information
<b>Tag</b>	C-His
<b>Accession #</b>	Q12874
<b>Host</b>	E.coli
<b>Species</b>	Human
<b>Molecular Mass</b>	44 kDa
<b>Buffer</b>	Phosphate buffered saline
<b>Shipping</b>	The product is shipped on dry ice/polar packs. Upon receipt, store it immediately at the temperature listed below.
<b>Stability&amp;Storage</b>	Store at $\leq -70^{\circ}\text{C}$ , stable for 6 months after receipt. Store at $\leq -70^{\circ}\text{C}$ , stable for 3 months under sterile conditions after opening. Please minimize freeze-thaw cycles.
<b>Sequence</b>	VMAKEMLTKK STLRDQINSD HRTRAMQDRY MEVSGNLRDL YDDKDGLRKE ELNAISGPNE FAEFYNRLKQ IKEFHRKHPN EICVPMSVEF EELLKAREN SEEAQNLVEF TDEEGYGRYL DLHDCYLKYI NLKASEKLDY ITYLSIFDQL FDIPKERKNA EYKRYLEMILL EYLQDYTDV KPLQDQNELF GKIQAEFK WENGTFPGWP KETSSALTHA GAHLDLSAFS SWEELASLGL DRLKSALLAL GLKCGGTL EE RAQRLFSTKG KSLES LDTSL FAKNP KSKGT KR DTERNKDI AFLEAQIY EY VEILGEQRHL THENVQRKQA RTGEER EEEEE EQISESESE

## image



## Background

<b>Alternative Names</b>	Splicing factor 3A subunit 3 (SF3a60) (Spliceosome-associated protein 61) (SAP 61)
--------------------------	--

**Product Name: Recombinant Human SF3A3 Protein**  
**Catalog #: PEH4892**

---



**Note**

For Research Use Only , Not for Diagnostic Use