


Product Name: Recombinant Human DNAJC3 Protein
Catalog #: PEH4894



Summary

Name	DNJC3_HUMAN
Purity	Greater than 85% as determined by reducing SDS-PAGE
Endotoxin level	Please contact with the lab for this information
Tag	C-His
Accession #	Q13217
Host	E.coli
Species	Human
Molecular Mass	44 kDa
Buffer	Phosphate buffered saline
Shipping	The product is shipped on dry ice/polar packs. Upon receipt, store it immediately at the temperature listed below.
Stability&Storage	Store at $\leq -70^{\circ}\text{C}$, stable for 6 months after receipt. Store at $\leq -70^{\circ}\text{C}$, stable for 3 months under sterile conditions after opening. Please minimize freeze-thaw cycles.
Sequence	RGHLLKQGK LDEAEDDFKK VLKSNPSENE EKEAQSQLIK SDEMQLRSQ ALNAFGSGDY TAAIAFLDKI LEVCVWDAEL RELRAECFIK EGEPKKAISD LKAASKLKND NTEAFYKIST LYYQLGDHEL SLSEVRECLK LDQDHKRCFA HYKQVKKLNK LIESAEELIR DGRYTDATSK YESVMKTEPS IAEYTVRSKE RICHCFKDE KPVEAIRVCS EVLQMEPDNV NALKDRAEAY LIEEMYDEAI QDYETAQEHN ENDQQIREGL EKAQRLKQS QKRDYYKILG VKRNAKKQEI IKAYRKLALQ WHPDNFQNEE EKKKAEEKFI DIAAAKEVLS DPEMRKKFDD GEDPLDAESQ

image



Background

Alternative Names	DnaJ homolog subfamily C member 3 (Endoplasmic reticulum DNA J domain-
--------------------------	--

Product Name: Recombinant Human DNAJC3 Protein
Catalog #: PEH4894



containing protein 6) (ER-resident protein ERdj6) (ERdj6) (Interferon-induced, double-stranded RNA-activated protein kinase inhibitor) (Protein kinase inhibitor of 58 kDa) (Protein kinase inhibitor p58)

Note

For Research Use Only , Not for Diagnostic Use