


Product Name: Recombinant Human SNW1 Protein
Catalog #: PEH4901



Summary

Name	SNW1_HUMAN
Purity	Greater than 85% as determined by reducing SDS-PAGE
Endotoxin level	Please contact with the lab for this information
Tag	C-His
Accession #	Q13573
Host	E.coli
Species	Human
Molecular Mass	16 kDa
Buffer	Phosphate buffered saline
Shipping	The product is shipped on dry ice/polar packs. Upon receipt, store it immediately at the temperature listed below.
Stability&Storage	Store at $\leq -70^{\circ}\text{C}$, stable for 6 months after receipt. Store at $\leq -70^{\circ}\text{C}$, stable for 3 months under sterile conditions after opening. Please minimize freeze-thaw cycles.
Sequence	IARQGQSKDK VIYSKYTDLV PKEVMNADDP DLQRPDEEAI KEITEKTRVA LEKSVSQKVA AAMPVRAADK LAPAQYIRYT PSQQGVAFNS GAKQRVIRMV EMQKDPMEPP RFKINKKIPR

image



Background

Alternative Names	SNW domain-containing protein 1 (Nuclear protein SkiP) (Nuclear receptor coactivator NCoA-62) (Ski-interacting protein)
--------------------------	---

Note

For Research Use Only , Not for Diagnostic Use