

Product Name: Recombinant Human DMC1 Protein
Catalog #: PEH4906



Summary

Name	DMC1_HUMAN
Purity	Greater than 85% as determined by reducing SDS-PAGE
Endotoxin level	Please contact with the lab for this information
Tag	C-His
Accession #	Q14565
Host	E.coli
Species	Human
Molecular Mass	38 kDa
Buffer	Lyophilized from Phosphate buffered saline
Shipping	The product is shipped on dry ice/polar packs. Upon receipt, store it immediately at the temperature listed below.
Stability&Storage	Store at $\leq -70^{\circ}\text{C}$, stable for 12 months after receipt. Store at $\leq -70^{\circ}\text{C}$, stable for 3 months under sterile conditions after opening. Please minimize freeze-thaw cycles.
Sequence	PGFQDEEESL FQDIDLLQKH GINVADIKKL KSVGICTIKG IQMTTRRALC NVKGLSEAKV DKIKEAANKL IEPGFLTAFE YSEKRKMVFH ITTGSQEFDK LLGGGIESMA ITEAFGEFRT GKTQLSHTLC VTAQLPGAGG YPGGKIIFID TENTFRPDRL RDIADRFNVD HDAVLDNVLY ARAYTSEHQM ELLDYVAAKF HEEAGIFKLL IIDSIMALFR VDFSGRGELA ERQQKLAQML SRLQKISEEY NVAVFVTNQM TADPGATMTF QADPKKPIGG HILAHASTTR ISLRKGRGEL RIAKIYDSPE MPENEATFAI

image



Background

Alternative Names Meiotic recombination protein DMC1/LIM15 homolog

Product Name: Recombinant Human DMC1 Protein
Catalog #: PEH4906



Note

For Research Use Only , Not for Diagnostic Use