


Product Name: Recombinant Human ZNF235 Protein
Catalog #: PEH4908



Summary

Name	ZN235_HUMAN
Purity	Greater than 85% as determined by reducing SDS-PAGE
Endotoxin level	Please contact with the lab for this information
Tag	C-His
Accession #	Q14590
Host	E.coli
Species	Human
Molecular Mass	22 kDa
Buffer	Lyophilized from Phosphate buffered saline
Shipping	The product is shipped on dry ice/polar packs. Upon receipt, store it immediately at the temperature listed below.
Stability&Storage	Store at $\leq -70^{\circ}\text{C}$, stable for 12 months after receipt. Store at $\leq -70^{\circ}\text{C}$, stable for 3 months under sterile conditions after opening. Please minimize freeze-thaw cycles.
Sequence	SQFPKHHDSP CQVGAGESIQ ASVDDNCLVN HIGDHSSIIH NQEFPTGKVP NSWSKIYLN TQNYQRSCKQ TQMKNKLCIF APYVDIFSCI SHHHDDNIVH KRDKVHSNSD CGKDTLKVSP LTQRSIHTGQ KTYQGNECEE AFNDSSSLEL HKQVHLGKKS PACSTHEKDT SYSSGIPVQQ

image



Background

Alternative Names Zinc finger protein 235 (Zinc finger protein 270) (Zinc finger protein 93 homolog) (Zfp-93) (Zinc finger protein HZF6)

Note

For Research Use Only , Not for Diagnostic Use

Product Name: Recombinant Human ZNF235 Protein
Catalog #: PEH4908

