


Product Name: Recombinant Human SMC1A Protein
Catalog #: PEH4911



Summary

Name	SMC1A_HUMAN
Purity	Greater than 85% as determined by reducing SDS-PAGE
Endotoxin level	Please contact with the lab for this information
Tag	C-His
Accession #	Q14683
Host	E.coli
Species	Human
Molecular Mass	56 kDa
Buffer	Lyophilized from Phosphate buffered saline
Shipping	The product is shipped on dry ice/polar packs. Upon receipt, store it immediately at the temperature listed below.
Stability&Storage	Store at $\leq -70^{\circ}\text{C}$, stable for 12 months after receipt. Store at $\leq -70^{\circ}\text{C}$, stable for 3 months under sterile conditions after opening. Please minimize freeze-thaw cycles.
Sequence	ELASKNKEIE KDKKRMDKVE DELKEKKEL GKMMREQQIQI EKEIKEKDSE LNQKRPQYIK AKENTSHKIK KLEAAKSLQ NAQKHYKKRK GDMDELEKEM LSVEKARQEF EERMEEESQS QGRDLTLEEN QVKKYHRLKE EASKRAATLA QELEKFNRDQ KADQDRDLLE ERKKVETEAQ IKQKLREIEE NQKRIEKL YITTSKQSLE EQKKLEGELT EEVEMAKRRI DEINKELNQV MEQLGDARID RQESSRQQRK AEIMESIKRL YPGSVYGRIL DLCQPTQKKY QIAVTKVLGK NMDAIIVDSE KTGRDCIQYI KEQRGEPETF LPLDYLEVKP TDEKLRELKG AKLVIDVIRY EPPHIKKALQ YACGNALVCD NVEDARRIAF GGHQRHKTVA LDGTLFQKSG VISGGASDLK AKARRWDEKA VDKLKEKKER LTELKEQMK

image



Background

Product Name: Recombinant Human SMC1A Protein
Catalog #: PEH4911



Alternative Names

Structural maintenance of chromosomes protein 1A (SMC protein 1A)
(SMC-1-alpha) (SMC-1A) (Sb1.8)

Note

For Research Use Only , Not for Diagnostic Use