


Product Name: Recombinant Human NCOA6 Protein
Catalog #: PEH4912



Summary

Name	NCOA6_HUMAN
Purity	Greater than 85% as determined by reducing SDS-PAGE
Endotoxin level	Please contact with the lab for this information
Tag	C-His
Accession #	Q14686
Host	E.coli
Species	Human
Molecular Mass	16 kDa
Buffer	Phosphate buffered saline
Shipping	The product is shipped on dry ice/polar packs. Upon receipt, store it immediately at the temperature listed below.
Stability&Storage	Store at $\leq -70^{\circ}\text{C}$, stable for 6 months after receipt. Store at $\leq -70^{\circ}\text{C}$, stable for 3 months under sterile conditions after opening. Please minimize freeze-thaw cycles.
Sequence	DIYTSLCSST MEDSEMDFDS GLEDDDTKSD SILEDSTIFV AFKGNIDDKD FKWKLDAILK NVPNLLHMES SKLKVQKVEP WNSVRVTFNI PREAAERLRI LAQSNNQQLR DLGILSVQIE GEGAINLALA

image



Background

Alternative Names	Nuclear receptor coactivator 6 (Activating signal cointegrator 2) (ASC-2) (Amplified in breast cancer protein 3) (Cancer-amplified transcriptional coactivator ASC-2) (Nuclear receptor coactivator RAP250) (NRC RAP250) (Nuclear receptor-activating protein, 250 kDa) (Peroxisome proliferator-activated receptor-interacting protein) (PPAR-interacting protein) (PRIP) (Thyroid hormone receptor-binding protein)
--------------------------	---

Product Name: Recombinant Human NCOA6 Protein
Catalog #: PEH4912



Note

For Research Use Only , Not for Diagnostic Use